



FountainCircle
Medical Suppliers

OxyPure

- 1-5 L/min
- Lightweight and Robust Design
- Suitable for Most Oxygen Prescriptions
- Fast and Easy Maintenance
- Oxygen and Performance Monitor



www.fcmedical.co.za



155 Anderson Street
Brooklyn - Pretoria 0181

OxyPure


OXYGEN CONCENTRATOR

Specification	Model: OxyPure
Oxygen Flow	1 to 5 Litres per minute
Delivery Method	Continuous
Sound Levels (dBA)	40 dBA
Weight (kg)	16 Kgs
Size (LxWxH)	36x29x58 cm
Cord Length (m)	2 m
Electrical Requirements	115V/60Hz or 220V/50Hz or 230V/60Hz

Up to 5 litre Flow



Humidifier



Easy to read display



Removable Filter

